

Fish Caspase 3, CASP3 ELISA Kit | CASP3 elisa kit

Fish Caspase 3, CASP3 ELISA Kit

CATALOG #
AAA123057

REACTIVITY
Fish


PRODUCT SYNONYMS

Caspase 3; CASP3; Fish Caspase 3; CASP3 ELISA Kit; CASP3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Samples	Serum, Plasma, Cell Culture Supernates, Cell Lysates, Tissue Homogenates
Assay Type	Quantitative Sandwich
Detection Range	10ng/L-2000ng/L
Sensitivity	5.01ng/L
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CASP3 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Fish Cardiac Troponin 1 (also known as cTn-1) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish cTn-1 antibody. cTn-1 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish cTn-1 Antibody is added and binds to cTn-1 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated cTn-1 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish cTn-1. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fish-elisa-kits/caspase-3-casp3/123057>