

Human Dimethylaniline monooxygenase [N-oxide-forming] 3, FMO3 ELISA Kit | FMO3 elisa kit

Human Dimethylaniline monooxygenase [N-oxide-forming] 3, FMO3 ELISA Kit

CATALOG # AAA123507	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Dimethylaniline monooxygenase [N-oxide-forming] 3; FMO3; Human Dimethylaniline monooxygenase [N-oxide-forming] 3; FMO3 ELISA Kit; FMO3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	571
Samples	Serum, plasma, cell culture supernates, ascites, tissue homogenates or other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	10-2000ng/L
Sensitivity	4.36ng/L
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FMO3 elisa kit	Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Dimethylaniline monooxygenase [N-oxide-forming] 3 (also known as FMO3) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids. Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human FMO3 antibody. FMO3 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human FMO3 Antibody is added and binds to FMO3 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated FMO3 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human FMO3. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	1241854429
NCBI Official Full Name	dimethylaniline monooxygenase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/dimethylaniline-monooxygenase-n-oxide-forming-3-fmo3/123507>