

Human Cluster of Differentiation 57, CD57 ELISA Kit | CD57 elisa kit

Human Cluster of Differentiation 57, CD57 ELISA Kit

CATALOG #
AAA123829

REACTIVITY
Human


PRODUCT SYNONYMS

Cluster of Differentiation 57; CD57; Human Cluster of Differentiation 57; CD57 ELISA Kit; CD57 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	0.1-40ng/ml
Sensitivity	0.044ng/ml
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store the reagents at 2-8 degree C. For long term storage refer to the expiration date keep it at -20 degree C. Avoid repeated thaw cycles. If individual reagents are opened it is recommended that the kit be used within 1 month.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD57 elisa kit	Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human CD57 antibody. CD57 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human CD57 Antibody is added and binds to CD57 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated CD57 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human CD57. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm. Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Cluster of Differentiation 57 (also known as CD57) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cluster-of-differentiation-57-cd57/123829>