

Human Proto-Oncogene Mas, MAS1 ELISA Kit | MAS1 elisa kit

Human Proto-Oncogene Mas, MAS1 ELISA Kit

CATALOG #
AAA123897

REACTIVITY
Human


PRODUCT SYNONYMS

Proto-Oncogene Mas; MAS1; Human Proto-Oncogene Mas; MAS1 ELISA Kit; MAS1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	15-3000ng/L
Sensitivity	6.31ng/L
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAS1 elisa kit

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human MAS1 antibody. MAS1 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human MAS1 Antibody is added and binds to MAS1 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated MAS1 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human MAS1. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Proto-Oncogene Mas (also known as MAS1) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/proto-oncogene-mas-mas1/123897>