

Human Trace amine-associated receptor 1, TAAR1 ELISA Kit | TAAR1 elisa kit

Human Trace amine-associated receptor 1, TAAR1 ELISA Kit

CATALOG #
AAA123916

REACTIVITY
Human


PRODUCT SYNONYMS

Trace amine-associated receptor 1; TAAR1; Human Trace amine-associated receptor 1; TAAR1 ELISA Kit; TAAR1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	0.19-12ng/ml
Sensitivity	0.1ng/ml
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TAAR1 elisa kit

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human TAAR1 antibody. TAAR1 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human TAAR1 Antibody is added and binds to TAAR1 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated TAAR1 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human TAAR1. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Trace amine-associated receptor 1 (also known as TAAR1) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/trace-amine-associated-receptor-1-taar1/123916>