

Mouse inducible nitric oxide synthase (iNOS) ELISA Kit | iNOS elisa kit

Mouse inducible nitric oxide synthase (iNOS) ELISA Kit

CATALOG # AAA12478	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

inducible nitric oxide synthase (iNOS); Mouse inducible nitric oxide synthase (iNOS) ELISA Kit; iNOS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	1147
Assay Type	Sandwich
Detection Range	100 U/mL-1.56 U/mL
Sensitivity	Up to 0.5 U/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for iNOS elisa kit	Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Mouse iNOS monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	2935553
Molecular Weight	130,628 Da
NCBI Official Full Name	inducible nitric oxide synthase

UniProt Protein Name	Nitric oxide synthase, inducible
UniProt Gene Name	Nos2
UniProt Synonym Gene Names	Inducible NOS; iNOS
UniProt Entry Name	NOS2_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/inducible-nitric-oxide-synthase-inos/12478>