

Mouse 5-Hydroxyindoleacetic acid (5-HIAA) ELISA Kit | 5-HIAA elisa kit

Mouse 5-Hydroxyindoleacetic acid (5-HIAA) ELISA Kit

CATALOG # AAA12499	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

5-Hydroxyindoleacetic acid (5-HIAA); Mouse 5-Hydroxyindoleacetic acid (5-HIAA) ELISA Kit; 5-HIAA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	580
Assay Type	Sandwich
Detection Range	200 ng/mL-3.12 ng/mL
Sensitivity	Up to 0.6 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for 5-HIAA elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Mouse 5-HIAA monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	299779097
NCBI GeneID	9384451
NCBI Accession #	YP_003734291.1

NCBI GenBank Nucleotide #	NC_014260.1
Molecular Weight	63,941 Da
NCBI Official Full Name	5 gene product
NCBI Official Symbol	5
NCBI Protein Information	gp5 baseplate hub subunit and tail lysozyme
UniProt Protein Name	Lysozyme
UniProt Gene Name	5
UniProt Entry Name	D7RMF4_9CAUD

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/5-hydroxyindoleacetic-acid-5-hiaa/12499>