

Mouse RANK/CD265 ELISA Kit | RANK/CD265 elisa kit

Mouse RANK/CD265 PicoKine Fast ELISA Kit

CATALOG # AAA125346	REACTIVITY Mouse
-------------------------------	----------------------------


PRODUCT SYNONYMS

RANK/CD265; Mouse RANK/CD265 PicoKine Fast ELISA Kit; Mouse RANK; Tumor necrosis factor receptor superfamily member 11A; Osteoclast differentiation factor receptor; ODFR; Receptor activator of NF-KB; CD265; Tnfrsf11a; Rank; RANK/CD265 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	Natural and recombinant Mouse Tnfrsf11a
Sequence Length	625
Samples	Cell culture supernates and serum
Detection Range	62.5pg/ml-4000pg/ml
Sensitivity	<2pg/ml
Intra-assay Precision	Intra-Assay Precision (Precision within an assay): Three samples of known concentration were tested on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-Assay Precision (Precision accross assays):Three samples of known concentration were tested in separate assays to assess inter-assay precision.
Preparation and Storage	Store at 4 degree C for 6 months. Store at -20 degree C for 12 months. Avoid multiple freeze-thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for RANK/CD265 elisa kit	Principle of the Assay: The Picokine™ Mouse Tnfrsf11a Pre-Coated ELISA (Enzyme-Linked Immunosorbent Assay) kit is a solid phase immunoassay specially designed to measure Mouse Tnfrsf11a with a 96-well strip plate that is pre-coated with antibody specific for Tnfrsf11a. The detection antibody is a biotinylated antibody specific for Tnfrsf11a. The capture antibody is monoclonal antibody from rat, the detection antibody is polyclonal antibody from goat. The kit contains recombinant Mouse Tnfrsf11a with immunogen: Expression system for standard: NSO; Immunogen sequence: Q30-P213. The kit is analytically validated with ready to use reagents. To measure Mouse Tnfrsf11a, add standards and samples to the wells, then add the biotinylated detection antibody. Wash the wells with PBS or TBS buffer, and add Avidin-Biotin-Peroxidase Complex (ABC-HRP). Wash away the unbounded ABC-HRP with PBS or TBS buffer and add TMB. TMB is substrate to HRP and will be catalyzed to produce a blue color product, which changes into yellow after adding acidic stop solution. The density of the yellow product is linearly propotional to Mouse Tnfrsf11a in the sample. Read the density of the yellow product in each well using a plate reader, and benchmark the sample wells' readings against the standard curve to determine the concentration of Mouse Tnfrsf11a in the sample.!!Background/Introduction: Receptor Activator of Nuclear Factor kappa B(RANK), also known as TRANCE Receptor, is a type I membrane protein that is expressed on the surface of osteoclasts and is involved in their activation upon ligand binding. RANK is a recently described TNF receptor family member, and its ligand, RANKL, promote survival of dendritic cells and differentiation of osteoclasts. RANK contains 383 amino acids in its intracellular domain(residues 234-616), which contain three putative TRAF-binding domains(termed I, II, and III). RANK interacts with various TRAFs through distinct motifs and activates NF-kappaB via a novel TRAF6 interaction motif, which then activates NIK, thus leading to NF-kappaB activation, whereas RANK most likely activates JNK through a TRAF2-interacting region in RANK.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	17380323
NCBI GeneID	21934
UniProt Accession #	O35305
NCBI Official Full Name	Tumor necrosis factor receptor superfamily member 11A
NCBI Official Synonym Full Names	tumor necrosis factor receptor superfamily, member 11a, NFkB activator
NCBI Official Symbol	Tnfrsf11a
NCBI Official Synonym Symbols	OFE; ODFR; Rank; Ly109; mRANK; TRANCE-R
NCBI Protein Information	tumor necrosis factor receptor superfamily member 11A
UniProt Protein Name	Tumor necrosis factor receptor superfamily member 11A
UniProt Gene Name	Tnfrsf11a
UniProt Synonym Gene Names	Rank; ODFR
UniProt Entry Name	TNR11_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.