

Human Adenosine Deaminase (ADA) ELISA Kit | ADA elisa kit

Human Adenosine Deaminase (ADA) ELISA Kit

CATALOG # AAA12540	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Adenosine Deaminase (ADA); Human Adenosine Deaminase (ADA) ELISA Kit; ADA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	363
Samples	Human serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	Up to 0.5ng/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADA elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human ADA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	47078295
NCBI GeneID	100

NCBI Accession #	NP_000013.2
NCBI GenBank Nucleotide #	NM_000022.2
Molecular Weight	40,764 Da
NCBI Official Full Name	adenosine deaminase
NCBI Official Synonym Full Names	adenosine deaminase
NCBI Official Symbol	ADA
NCBI Protein Information	adenosine deaminase; adenosine aminohydrolase
UniProt Protein Name	Adenosine deaminase
UniProt Gene Name	ADA
UniProt Synonym Gene Names	ADA1
UniProt Entry Name	ADA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/adenosine-deaminase-ada/12540>