

Rat 50% complement hemolysis (CH50) ELISA Kit | CH50 elisa kit

Rat 50% complement hemolysis (CH50) ELISA Kit

CATALOG #
AAA12551

REACTIVITY
Rat


PRODUCT SYNONYMS

50% complement hemolysis (CH50); Rat 50% complement hemolysis (CH50) ELISA Kit; CH50 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	1741
Assay Type	Sandwich
Detection Range	100 U/mL-1.56 U/mL
Sensitivity	Up to 0.5 U/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CH50 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Rat CH50 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	262050656
NCBI GeneID	280678

NCBI Accession #	NP_001159957.1
NCBI GenBank Nucleotide #	NM_001166485.1
Molecular Weight	101,885 Da
NCBI Official Full Name	complement C4
NCBI Official Symbol	C4A
NCBI Official Synonym Symbols	C4
NCBI Protein Information	complement C4
UniProt Protein Name	Complement C4
UniProt Gene Name	C4
UniProt Entry Name	CO4_BOVIN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/50-complement-hemolysis-ch50/12551>