

Human Matrix Gla Protein (MGP) ELISA Kit | MGP elisa kit

Human Matrix Gla Protein (MGP) ELISA Kit

CATALOG #
AAA12560

REACTIVITY
Human


PRODUCT SYNONYMS

Matrix Gla Protein (MGP); Human Matrix Gla Protein (MGP) ELISA Kit; MGP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	103
Assay Type	Sandwich
Detection Range	1000 pg/mL-15.6 pg/mL
Sensitivity	Up to 5 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MGP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human MGP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

[187591](#)

NCBI GeneID	4256
Molecular Weight	15,325 Da
NCBI Official Full Name	matrix Gla protein
NCBI Official Synonym Full Names	matrix Gla protein
NCBI Official Symbol	MGP
NCBI Official Synonym Symbols	NTI; GIG36; MGLAP
NCBI Protein Information	matrix Gla protein; cell growth-inhibiting gene 36 protein
UniProt Protein Name	Matrix Gla protein
UniProt Gene Name	MGP
UniProt Synonym Gene Names	MGLAP; GIG36; MGP
UniProt Entry Name	MGP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/matrix-gla-protein-mgp/12560>