

Human Glial Fibrillary Acidic Protein (GFAP) ELISA Kit | GFAP elisa kit

Human Glial Fibrillary Acidic Protein (GFAP) ELISA Kit

CATALOG # AAA12611	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Glial Fibrillary Acidic Protein (GFAP); Human Glial Fibrillary Acidic Protein (GFAP) ELISA Kit; GFAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	432
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GFAP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human GFAP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	251802
NCBI GeneID	2670

Molecular Weight	49,508 Da
NCBI Official Full Name	glial fibrillary acidic protein
NCBI Official Synonym Full Names	glial fibrillary acidic protein
NCBI Official Symbol	GFAP
NCBI Protein Information	glial fibrillary acidic protein
UniProt Protein Name	Glial fibrillary acidic protein
UniProt Gene Name	GFAP
UniProt Synonym Gene Names	GFAP
UniProt Entry Name	GFAP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/glial-fibrillary-acidic-protein-gfap/12611>