

Human IGFBP7/Igfbp Rp1 ELISA Kit | IGFBP7 elisa kit

Human IGFBP7/Igfbp Rp1 PicoKine Quick ELISA Kit

CATALOG # AAA126825	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

IGFBP7/Igfbp Rp1; Human IGFBP7/Igfbp Rp1 PicoKine Quick ELISA Kit; IGFBP7 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	Natural and recombinant Human IGFBP7
Assay Type	Sandwich
Samples	Cell culture supernatants, serum, plasma (heparin, EDTA) and urine
Detection Range	625 pg/ml - 40,000 pg/ml
Sensitivity	<300 pg/ml
Applicable Applications for IGFBP7 elisa kit	ELISA
Intra-assay Precision	Intra-Assay Precision (Precision within an assay): Three samples of known concentration were tested on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-Assay Precision (Precision across assays): Three samples of known concentration were tested in separate assays to assess inter-assay precision.
Preparation and Storage	Store at 4 degree C for 6 months, at -20 degree C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

IMAGES



Standard Curve (Sample)

Human IGFBP7 PicoKine ELISA Kit standard curve

ADDITIONAL INFORMATION

Related Product Information for IGFBP7 elisa kit	Background: Insulin-like growth factor-binding protein 7 is a protein that in humans is encoded by the IGFBP7 gene. The gene was mapped to chromosome 4q12. The major function of the protein is the regulation of availability of insulin-like growth factors(IGFs) in tissue as well as in modulating IGF binding to its receptors. IGFBP7 binds to IGF with high affinity. It also stimulates cell adhesion. The protein is implicated in some cancers. Principle of the Assay: The Quick Picokine™ Human IGFBP7 Pre-Coated ELISA (Enzyme-Linked Immunosorbent Assay) kit is a solid phase immunoassay specially designed to measure Human IGFBP7 in cell culture supernatants, serum, plasma (heparin, EDTA) and urine. It uses our proprietary Quick ELISA technology. Quick ELISA technology employs capture antibodies conjugated to an affinity tag that is recognized by the polyclonal antibody used to coat our Quick ELISA plates. This approach to sandwich ELISA allows the formation of the antibody-analyte sandwich complex in a single step, significantly reducing assay time. The kit includes Human IGFBP7 protein as standards. To measure Human IGFBP7, add standards and samples to the wells, then add antibody cocktail. Wash the wells with PBS or TBS buffer, and add TMB. TMB is an HRP substrate and will be catalyzed to produce a blue color product, which changes into yellow after adding the acidic stop solution. The absorbance of the yellow product is linearly proportional to Human IGFBP7 in the sample. Read the absorbance of the yellow product in each well using a plate reader, and benchmark the sample wells' readings against the standard curve to determine the concentration of Human IGFBP7 in the sample.
Product Categories/Family for IGFBP7 elisa kit	PicoKine ELISA Kits, Fast ELISA Kits

NCBI AND UNIPROT INFORMATION

NCBI GI #	23396609
NCBI GeneID	3490
UniProt Accession #	Q16270; Q07822; Q53YE6; Q9UCA8; B4E1N2; B7Z9W7
Molecular Weight	282
NCBI Official Full Name	Insulin-like growth factor-binding protein 7
NCBI Official Synonym Full Names	insulin-like growth factor binding protein 7
NCBI Official Symbol	IGFBP7
NCBI Official Synonym Symbols	AGM; PSF; TAF; FSTL2; IBP-7; MAC25; IGFBP-7; RAMSVPS; IGFBP-7v; IGFBPRP1
NCBI Protein Information	insulin-like growth factor-binding protein 7; IGFBP-rP1; angiomodulin; IGF-binding protein 7; PGI2-stimulating factor; tumor-derived adhesion factor; prostacyclin-stimulating factor
UniProt Protein Name	Insulin-like growth factor-binding protein 7
UniProt Gene Name	IGFBP7
UniProt Synonym Gene Names	MAC25; PSF; IBP-7; IGF-binding protein 7; IGFBP-7; TAF
UniProt Entry Name	IBP7_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.