

Human Dihydrotestosterone (DHT) ELISA Kit | DHT elisa kit

Human Dihydrotestosterone (DHT) ELISA Kit

CATALOG #
AAA12907

REACTIVITY
Human


PRODUCT SYNONYMS

Dihydrotestosterone (DHT); Human Dihydrotestosterone (DHT) ELISA Kit; DHT elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	381
Assay Type	Sandwich
Detection Range	2000 pg/mL-31.2 pg/mL
Sensitivity	Up to 12 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DHT elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human DHT monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

1903379

Molecular Weight	43,009 Da
NCBI Official Full Name	differentially repressed by testosterone and dihydrotestosterone
UniProt Protein Name	Protein NDRG1
UniProt Gene Name	Ndr g1
UniProt Synonym Gene Names	Ndr1; NdrI; Tdd5; Protein Ndr1
UniProt Entry Name	NDRG1_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/dihydrotestosterone-dht/12907>