

Human Bovine Serum Albumin (BSA) ELISA Kit | BSA elisa kit

Human Bovine Serum Albumin (BSA) ELISA Kit

CATALOG #
AAA12913

REACTIVITY
Human


PRODUCT SYNONYMS

Bovine Serum Albumin (BSA); Human Bovine Serum Albumin (BSA) ELISA Kit; BSA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	607
Assay Type	Sandwich
Detection Range	200 ug/mL-3.12 ug/mL
Sensitivity	Up to 0.6 ug/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BSA elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human BSA monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	3336842
NCBI GeneID	280717

UniProt Accession #	P02769; O02787; P04277; Q3SZR2; A5PJX3
Molecular Weight	69,293 Da
NCBI Official Full Name	bovine serum albumin
NCBI Official Synonym Full Names	albumin
NCBI Official Symbol	ALB
NCBI Protein Information	serum albumin; BSA
UniProt Protein Name	Serum albumin
UniProt Gene Name	ALB
UniProt Entry Name	ALBU_BOVIN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/bovine-serum-albumin-bsa/12913>