

Human Ubiquitin Conjugating Enzyme E2I (UBE2I) ELISA Kit | UBE2I elisa kit

Human Ubiquitin Conjugating Enzyme E2I (UBE2I) ELISA Kit

CATALOG #
AAA129529

REACTIVITY
Human


PRODUCT SYNONYMS

Ubiquitin Conjugating Enzyme E2I (UBE2I); Human Ubiquitin Conjugating Enzyme E2I (UBE2I) ELISA Kit; P18; UBC9; UBCE9; SUMO-protein ligase; Ubiquitin carrier protein 9; Ubiquitin carrier protein I; Ubiquitin-protein ligase I; SUMO-conjugating enzyme UBC9; UBE2I elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Ubiquitin Conjugating Enzyme E2I (UBE2I). No significant cross-reactivity or interference between Ubiquitin Conjugating Enzyme E2I (UBE2I) and analogues was observed.
Sequence Length	158
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	0.156-10ng/mL
Sensitivity	0.057ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Ubiquitin Conjugating Enzyme E2I (UBE2I) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Ubiquitin Conjugating Enzyme E2I (UBE2I) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UBE2I elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Ubiquitin Conjugating Enzyme E2I (UBE2I). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Ubiquitin Conjugating Enzyme E2I (UBE2I). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Ubiquitin Conjugating Enzyme E2I (UBE2I), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Ubiquitin Conjugating Enzyme E2I (UBE2I) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for UBE2I elisa kit	Enzyme & Kinase

NCBI AND UNIPROT INFORMATION

NCBI GI #	532691743
NCBI GeneID	101955823
NCBI Accession #	NP_001269197.1
NCBI GenBank Nucleotide #	NM_001282268.1
UniProt Accession #	Q2EF73
Molecular Weight	18,007 Da
NCBI Official Full Name	SUMO-conjugating enzyme UBC9
NCBI Official Symbol	Ube2i
NCBI Official Synonym Symbols	Ubc9
NCBI Protein Information	SUMO-conjugating enzyme UBC9
UniProt Protein Name	SUMO-conjugating enzyme UBC9
UniProt Gene Name	UBE2I
UniProt Synonym Gene Names	UBC9
UniProt Entry Name	UBC9_SPETR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/ubiquitin-conjugating-enzyme-e2i-ube2i/129529>