

Human Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) ELISA Kit | ALPPL2 elisa kit

Human Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) ELISA Kit

CATALOG #
AAA129721

REACTIVITY
Human


PRODUCT SYNONYMS

Alkaline Phosphatase; Placental Like Protein 2 (ALPPL2); Human Alkaline Phosphatase; Placental Like Protein 2 (ALPPL2) ELISA Kit; ALPG; GCAP; ALP-1; Alkaline phosphatase Nagao isozyme; Germ cell alkaline phosphatase; Placental alkaline phosphatase-like; ALPPL2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2). No significant cross-reactivity or interference between Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) and analogues was observed.
Sequence Length	532
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.117ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ALPPL2 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for ALPPL2 elisa kit	Enzyme & Kinase

NCBI AND UNIPROT INFORMATION

NCBI GI #	157266296
NCBI GeneID	251
NCBI Accession #	NP_112603.2
NCBI GenBank Nucleotide #	NM_031313.2
Molecular Weight	57,377 Da
NCBI Official Full Name	alkaline phosphatase, placental-like preproprotein
NCBI Official Synonym Full Names	alkaline phosphatase, placental-like 2
NCBI Official Symbol	ALPPL2
NCBI Official Synonym Symbols	ALPG; GCAP; ALPPL
NCBI Protein Information	alkaline phosphatase, placental-like
UniProt Protein Name	Alkaline phosphatase, placental-like
UniProt Gene Name	ALPPL2
UniProt Synonym Gene Names	ALPPL; GCAP; PLAP-like
UniProt Entry Name	PPBN_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.