

Porcine Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) ELISA Kit | Vaspin elisa kit

Porcine Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) ELISA Kit

CATALOG #
AAA12985

REACTIVITY
Porcine


PRODUCT SYNONYMS

Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin); Porcine Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) ELISA Kit; Vaspin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	No cross-reaction with other factors.
Sequence Length	414
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Vaspin elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Porcine Vaspin monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	32816048
Molecular Weight	47,175 Da
NCBI Official Full Name	visceral adipose-specific SERPIN
UniProt Protein Name	Serpin A12
UniProt Gene Name	SERPINA12
UniProt Synonym Gene Names	Vaspin
UniProt Entry Name	SPA12_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/visceral-adipose-tissue-derived-serine-protease-inhibitor-vaspin/12985>