

Human Filaggrin (FLG) CLIA Kit | FLG clia kit

Human Filaggrin (FLG) CLIA Kit

CATALOG #
AAA129943

REACTIVITY
Human


PRODUCT SYNONYMS

Filaggrin (FLG); Human Filaggrin (FLG) CLIA Kit; Profilaggrin; FLG clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Filaggrin (FLG). No significant cross-reactivity or interference between Filaggrin (FLG) and analogues was observed.
Sequence Length	4061
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	27.4-20,000pg/mL
Sensitivity	9.5pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Filaggrin (FLG) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Filaggrin (FLG) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FLG clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Filaggrin (FLG). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Filaggrin (FLG). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Filaggrin (FLG) level in the sample or standard.;
Product Categories/Family for FLG clia kit	Metabolic pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	60097902
NCBI GeneID	2312
NCBI Accession #	NP_002007.1
NCBI GenBank Nucleotide #	NM_002016.1
UniProt Accession #	P20930; Q01720; Q5T583; Q9UC71
Molecular Weight	435,170 Da
NCBI Official Full Name	filaggrin
NCBI Official Synonym Full Names	filaggrin
NCBI Official Symbol	FLG
NCBI Official Synonym Symbols	ATOD2
NCBI Protein Information	filaggrin
UniProt Protein Name	Filaggrin
UniProt Gene Name	FLG
UniProt Entry Name	FILA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.