

Human Complement Component 8g (C8g) CLIA Kit | C8g clia kit

Human Complement Component 8g (C8g) CLIA Kit

CATALOG #
AAA130231

REACTIVITY
Human


PRODUCT SYNONYMS

Complement Component 8g (C8g); Human Complement Component 8g (C8g) CLIA Kit; C8g clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Complement Component 8g (C8g). No significant cross-reactivity or interference between Complement Component 8g (C8g) and analogues was observed.
Sequence Length	237
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma and other biological fluids
Detection Range	6.86-5,000pg/mL
Sensitivity	2.69pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Complement Component 8g (C8g) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Complement Component 8g (C8g) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for C8g clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Complement Component 8g (C8g). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Complement Component 8g (C8g). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Complement Component 8g (C8g) level in the sample or standard.;
Product Categories/Family for C8g clia kit	Metabolic pathway ; Immune molecule ; Hematology

NCBI AND UNIPROT INFORMATION

NCBI GI #	158519861
NCBI GeneID	517356
NCBI Accession #	NP_001103546.1
NCBI GenBank Nucleotide #	NM_001110076.2
NCBI Official Full Name	complement component C8 gamma chain
NCBI Official Symbol	C8G
NCBI Protein Information	complement component C8 gamma chain

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/complement-component-8g-c8g/130231>