

Human Semaphorin 3A (SEMA3A) CLIA Kit | SEMA3A clia kit

Human Semaphorin 3A (SEMA3A) CLIA Kit

CATALOG #
AAA130239

REACTIVITY
Human


PRODUCT SYNONYMS

Semaphorin 3A (SEMA3A); Human Semaphorin 3A (SEMA3A) CLIA Kit; SEMAL; Hsema-I; Hsema-III; SEMA1; SEMAD; SEMAIII; SemD; Coll-1; Sema Domain; Immunoglobulin Domain(Ig); Transmembrane Domain(TM)and Short Cytoplasmic Domain 3A; SEMA3A clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Semaphorin 3A (SEMA3A). No significant cross-reactivity or interference between Semaphorin 3A (SEMA3A) and analogues was observed.
Sequence Length	771
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	13.7-10,000pg/mL
Sensitivity	5.1pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Semaphorin 3A (SEMA3A) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Semaphorin 3A (SEMA3A) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SEMA3A clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Semaphorin 3A (SEMA3A). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Semaphorin 3A (SEMA3A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Semaphorin 3A (SEMA3A) level in the sample or standard.;
Product Categories/Family for SEMA3A clia kit	CD & Adhesion molecule

NCBI AND UNIPROT INFORMATION

NCBI GI #	5174673
NCBI GeneID	10371
NCBI Accession #	NP_006071.1
NCBI GenBank Nucleotide #	NM_006080.2
UniProt Accession #	Q14563
Molecular Weight	88,889 Da
NCBI Official Full Name	semaphorin-3A
NCBI Official Synonym Full Names	sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A
NCBI Official Symbol	SEMA3A
NCBI Official Synonym Symbols	HH16; SemD; COLL1; SEMA1; SEMAD; SEMAL; coll-1; Hsema-I; SEMAIII; Hsema-III
NCBI Protein Information	semaphorin-3A
UniProt Protein Name	Semaphorin-3A
UniProt Gene Name	SEMA3A
UniProt Synonym Gene Names	SEMAD; Sema III
UniProt Entry Name	SEM3A_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.