

Human Epstein Barr Virus Induced Protein 3 (EBI3) CLIA Kit | EBI3 clia kit

Human Epstein Barr Virus Induced Protein 3 (EBI3) CLIA Kit

CATALOG #
AAA130276

REACTIVITY
Human


PRODUCT SYNONYMS

Epstein Barr Virus Induced Protein 3 (EBI3); Human Epstein Barr Virus Induced Protein 3 (EBI3) CLIA Kit; IL27B; IL27-B; Interleukin 27B; Epstein-Barr virus-induced gene 3 protein; EBV-induced gene 3 protein; EBI3 clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Epstein Barr Virus Induced Protein 3 (EBI3). No significant cross-reactivity or interference between Epstein Barr Virus Induced Protein 3 (EBI3) and analogues was observed.
Sequence Length	245
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	2.74-2,000pg/mL
Sensitivity	1.17pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Epstein Barr Virus Induced Protein 3 (EBI3) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Epstein Barr Virus Induced Protein 3 (EBI3) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EBI3 clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Epstein Barr Virus Induced Protein 3 (EBI3). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Epstein Barr Virus Induced Protein 3 (EBI3). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Epstein Barr Virus Induced Protein 3 (EBI3) level in the sample or standard.;
Product Categories/Family for EBI3 clia kit	Infection immunity

NCBI AND UNIPROT INFORMATION

NCBI GI #	659667890
NCBI GeneID	103380170
NCBI Accession #	NP_001281166.1
NCBI GenBank Nucleotide #	NM_001294237.1
NCBI Official Full Name	Epstein-Barr virus induced 3
NCBI Official Symbol	ebi3
NCBI Protein Information	Epstein-Barr virus induced 3

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/epstein-barr-virus-induced-protein-3-ebi3/130276>