

Canine Interferon Gamma (IFNg) CLIA Kit | IFNg clia kit

Canine Interferon Gamma (IFNg) CLIA Kit

CATALOG #
AAA130290

REACTIVITY
Canine


PRODUCT SYNONYMS

Interferon Gamma (IFNg); Canine Interferon Gamma (IFNg) CLIA Kit; IFN-G; IFG; IFI; INFr; IFN; Immune Interferon; IFNg clia kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	This assay has high sensitivity and excellent specificity for detection of Interferon Gamma (IFNg). No significant cross-reactivity or interference between Interferon Gamma (IFNg) and analogues was observed.
Sequence Length	166
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma and other biological fluids
Detection Range	1.37-1,000pg/mL
Sensitivity	0.65pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Interferon Gamma (IFNg) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Interferon Gamma (IFNg) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IFNg clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Interferon Gamma (IFNg). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Interferon Gamma (IFNg). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Interferon Gamma (IFNg) level in the sample or standard.;
Product Categories/Family for IFNg clia kit	Cytokine; Infection immunity

NCBI AND UNIPROT INFORMATION

NCBI GI #	56786138
NCBI GeneID	3458
NCBI Accession #	NP_000610.2
NCBI GenBank Nucleotide #	NM_000619.2
UniProt Accession #	P01579; Q53ZV4; B5BU88
Molecular Weight	19,348 Da
NCBI Official Full Name	interferon gamma
NCBI Official Synonym Full Names	interferon, gamma
NCBI Official Symbol	IFNG
NCBI Official Synonym Symbols	IFG; IFI
NCBI Protein Information	interferon gamma
UniProt Protein Name	Interferon gamma
UniProt Gene Name	IFNG
UniProt Synonym Gene Names	IFN-gamma
UniProt Entry Name	IFNG_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-clia-kits/interferon-gamma-ifng/130290>