

Mouse N-Acetyl Beta-D-Glucosaminidase (NAGase) CLIA Kit | NAGase clia kit

Mouse N-Acetyl Beta-D-Glucosaminidase (NAGase) CLIA Kit

CATALOG #
AAA130497

REACTIVITY
Mouse


PRODUCT SYNONYMS

N-Acetyl Beta-D-Glucosaminidase (NAGase); Mouse N-Acetyl Beta-D-Glucosaminidase (NAGase) CLIA Kit; MGEA5; MEA5; NCOAT; Hexosaminidase C; N-Acetylhexosaminidase; Hyaluronidase; Meningioma-Expressed Antigen 5; Nuclear Cytoplasmic O-GlcNAcase And Acetyltransferase; NAGase clia kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of N-Acetyl Beta-D-Glucosaminidase (NAGase). No significant cross-reactivity or interference between N-Acetyl Beta-D-Glucosaminidase (NAGase) and analogues was observed.
Sequence Length	627
Assay Type	Double-antibody Sandwich
Samples	Urine
Detection Range	68.6-50,000pg/mL
Sensitivity	29.4pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level N-Acetyl Beta-D-Glucosaminidase (NAGase) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level N-Acetyl Beta-D-Glucosaminidase (NAGase) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NAGase clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to N-Acetyl Beta-D-Glucosaminidase (NAGase). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to N-Acetyl Beta-D-Glucosaminidase (NAGase). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the N-Acetyl Beta-D-Glucosaminidase (NAGase) level in the sample or standard.;
Product Categories/Family for NAGase clia kit	Enzyme & Kinase; Kidney biomarker

NCBI AND UNIPROT INFORMATION

NCBI GI #	41393807
Molecular Weight	72,252 Da
NCBI Official Full Name	N-acetyl-beta-D-glucosaminidase, partial
UniProt Protein Name	N-acetyl-beta-D-glucosaminidase
UniProt Entry Name	Q6ST21_STRGN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-clia-kits/n-acetyl-beta-d-glucosaminidase-nagase/130497>