

Human Toll Like Receptor 9 (TLR9) CLIA Kit | TLR9 clia kit

Human Toll Like Receptor 9 (TLR9) CLIA Kit

CATALOG #
AAA130565

REACTIVITY
Human


PRODUCT SYNONYMS

Toll Like Receptor 9 (TLR9); Human Toll Like Receptor 9 (TLR9) CLIA Kit; CD289; TLR9 clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Toll Like Receptor 9 (TLR9). No significant cross-reactivity or interference between Toll Like Receptor 9 (TLR9) and analogues was observed.
Sequence Length	1032
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	2.74-2,000pg/mL
Sensitivity	0.99pg/mL
Application	Chemiluminescent immunoassay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Toll Like Receptor 9 (TLR9) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Toll Like Receptor 9 (TLR9) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TLR9 clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Toll Like Receptor 9 (TLR9). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Toll Like Receptor 9 (TLR9). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Toll Like Receptor 9 (TLR9) level in the sample or standard.;
Product Categories/Family for TLR9 clia kit	Immune molecule

NCBI AND UNIPROT INFORMATION

NCBI GI #	157057166
NCBI GeneID	81897
NCBI Accession #	NP_112455.2
NCBI GenBank Nucleotide #	NM_031178.2
UniProt Accession #	Q9EQU3; Q4L0K3; Q4L0K4; Q99MF2; Q99MQ8; F8VPN5
Molecular Weight	116,412 Da
NCBI Official Full Name	toll-like receptor 9
NCBI Official Synonym Full Names	toll-like receptor 9
NCBI Official Symbol	Tlr9
NCBI Protein Information	toll-like receptor 9
UniProt Protein Name	Toll-like receptor 9
UniProt Gene Name	Tlr9
UniProt Entry Name	TLR9_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/toll-like-receptor-9-tlr9/130565>