

Human Protein Kinase, AMP Activated Beta 1 (PRKAb1) ELISA Kit | PRKAb1 elisa kit

Human Protein Kinase, AMP Activated Beta 1 (PRKAb1) ELISA Kit

CATALOG #
AAA135763

REACTIVITY
Human


PRODUCT SYNONYMS

Protein Kinase; AMP Activated Beta 1 (PRKAb1); Human Protein Kinase; AMP Activated Beta 1 (PRKAb1) ELISA Kit; AMPK; HAMPKb; 5'-AMP-Activated Protein Kinase Subunit Beta-1 Non-Catalytic Subunit; Adenosine Monophosphate Activated Protein Kinase(AMPK); PRKAb1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Protein Kinase, AMP Activated Beta 1 (PRKAb1). No significant cross-reactivity or interference between Protein Kinase, AMP Activated Beta 1 (PRKAb1) and analogues was observed.
Sequence Length	260
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.124ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Protein Kinase, AMP Activated Beta 1 (PRKAb1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Protein Kinase, AMP Activated Beta 1 (PRKAb1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PRKAb1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Protein Kinase, AMP Activated Beta 1 (PRKAb1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Protein Kinase, AMP Activated Beta 1 (PRKAb1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Protein Kinase, AMP Activated Beta 1 (PRKAb1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Protein Kinase, AMP Activated Beta 1 (PRKAb1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for PRKAb1 elisa kit	Enzyme & Kinase; Metabolic pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	47550723
NCBI GeneID	335232
NCBI Accession #	NP_999878.1
NCBI GenBank Nucleotide #	NM_214713.1
NCBI Official Full Name	protein kinase, AMP-activated, beta 1 non-catalytic subunit, b
NCBI Official Synonym Full Names	protein kinase, AMP-activated, beta 1 non-catalytic subunit, b
NCBI Official Symbol	prkab1b
NCBI Official Synonym Symbols	prkab1; zgc:56652; zgc:76975; wu:fk93d05; wu:fw87e09
NCBI Protein Information	protein kinase, AMP-activated, beta 1 non-catalytic subunit, b

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.