

Human Vasodilator Stimulated Phosphoprotein (VASP) ELISA Kit | VASP elisa kit

Human Vasodilator Stimulated Phosphoprotein (VASP) ELISA Kit

CATALOG #
AAA136331

REACTIVITY
Human


PRODUCT SYNONYMS

Vasodilator Stimulated Phosphoprotein (VASP); Human Vasodilator Stimulated Phosphoprotein (VASP) ELISA Kit; VASP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Vasodilator Stimulated Phosphoprotein (VASP). No significant cross-reactivity or interference between Vasodilator Stimulated Phosphoprotein (VASP) and analogues was observed.
Sequence Length	380
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	0.156-10ng/mL
Sensitivity	0.059ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Vasodilator Stimulated Phosphoprotein (VASP) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Vasodilator Stimulated Phosphoprotein (VASP) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VASP elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Vasodilator Stimulated Phosphoprotein (VASP). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Vasodilator Stimulated Phosphoprotein (VASP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Vasodilator Stimulated Phosphoprotein (VASP), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Vasodilator Stimulated Phosphoprotein (VASP) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for VASP elisa kit	Signal Transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507869
NCBI GeneID	7408
NCBI Accession #	NP_003361.1
NCBI GenBank Nucleotide #	NM_003370.3
UniProt Accession #	P50552; Q6PIZ1; Q93035; B2RBT9
Molecular Weight	39,830 Da
NCBI Official Full Name	vasodilator-stimulated phosphoprotein
NCBI Official Synonym Full Names	vasodilator-stimulated phosphoprotein
NCBI Official Symbol	VASP
NCBI Protein Information	vasodilator-stimulated phosphoprotein
UniProt Protein Name	Vasodilator-stimulated phosphoprotein
UniProt Gene Name	VASP
UniProt Synonym Gene Names	VASP
UniProt Entry Name	VASP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/vasodilator-stimulated-phosphoprotein-vasp/136331>