

Rat Prostaglandin E Synthase, Microsomal (PTGES) ELISA Kit | PTGES elisa kit

Rat Prostaglandin E Synthase, Microsomal (PTGES) ELISA Kit

CATALOG #
AAA136692

REACTIVITY
Rat


PRODUCT SYNONYMS

Prostaglandin E Synthase; Microsomal (PTGES); Rat Prostaglandin E Synthase; Microsomal (PTGES) ELISA Kit; MGST1L1; MGST-IV; MGST1-L1; PGES; PIG12; PP102; PP1294; TP53I12; PGE Synthase; Tumor Protein p53 Inducible Protein 12; Glutathione S-Transferase 1-Like 1; MGST1-Like 1; PTGES elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	This assay has high sensitivity and excellent specificity for detection of Prostaglandin E Synthase, Microsomal (PTGES). No significant cross-reactivity or interference between Prostaglandin E Synthase, Microsomal (PTGES) and analogues was observed.
Sequence Length	152
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.115ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Prostaglandin E Synthase, Microsomal (PTGES) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Prostaglandin E Synthase, Microsomal (PTGES) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PTGES elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Prostaglandin E Synthase, Microsomal (PTGES). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Prostaglandin E Synthase, Microsomal (PTGES). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Prostaglandin E Synthase, Microsomal (PTGES), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Prostaglandin E Synthase, Microsomal (PTGES) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for PTGES elisa kit	Enzyme & Kinase; Metabolic pathway; Infection immunity

NCBI AND UNIPROT INFORMATION

NCBI GI #	4758910
NCBI GeneID	9536
NCBI Accession #	NP_004869.1
NCBI GenBank Nucleotide #	NM_004878.4
UniProt Accession #	O14684; O14900; Q5SZC0
Molecular Weight	17,102 Da
NCBI Official Full Name	prostaglandin E synthase
NCBI Official Synonym Full Names	prostaglandin E synthase
NCBI Official Symbol	PTGES
NCBI Official Synonym Symbols	PGES; MPGES; PIG12; PP102; PP1294; MGST-IV; MGST1L1; TP53I12; mPGES-1; MGST1-L1
NCBI Protein Information	prostaglandin E synthase
UniProt Protein Name	Prostaglandin E synthase
UniProt Gene Name	PTGES
UniProt Synonym Gene Names	MGST1L1; MPGES1; PGES; PIG12; MGST1-L1; MPGES-1
UniProt Entry Name	PTGES_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.