

## Human Trypsinogen Activation Peptide (TAP) ELISA Kit | TAP elisa kit

Human Trypsinogen Activation Peptide (TAP) ELISA Kit

CATALOG #  
**AAA137128**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Trypsinogen Activation Peptide (TAP); Human Trypsinogen Activation Peptide (TAP) ELISA Kit; TAP elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity             | This assay has high sensitivity and excellent specificity for detection of Trypsinogen Activation Peptide (TAP). No significant cross-reactivity or interference between Trypsinogen Activation Peptide (TAP) and analogues was observed.                                                                                                                                                                                                                                                  |
| Samples                 | Serum, plasma and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Assay Type              | Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Detection Range         | 0.25-16nmol/L The standard curve concentrations used for the ELISA's were 16nmol/L, 8nmol/L, 4nmol/L, 2nmol/L, 1nmol/L, 0.5nmol/L, 0.25nmol/L                                                                                                                                                                                                                                                                                                                                              |
| Sensitivity             | Typically less than 0.08nmol/L.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Organism Species        | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| ELISA Type              | Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Precision               | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Trypsinogen Activation Peptide (TAP) were tested 20 times on one plate, respectively.<br>Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Trypsinogen Activation Peptide (TAP) were tested on 3 different plates, 8 replicates in each plate.<br>$CV(\%) = SD/mean \times 100$<br>Intra-Assay: CV<10%<br>Inter-Assay: CV<12%                       |
| Preparation and Storage | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-----------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for TAP elisa kit | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Trypsinogen Activation Peptide (TAP). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Trypsinogen Activation Peptide (TAP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Trypsinogen Activation Peptide (TAP), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Trypsinogen Activation Peptide (TAP) in the samples is then determined by comparing the O.D. of the samples to the standard curve. |
|-----------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/trypsinogen-activation-peptide-tap/137128>