

Human Baculoviral IAP Repeat Containing Protein 6 (BIRC6) ELISA Kit | BIRC6 elisa kit

Human Baculoviral IAP Repeat Containing Protein 6 (BIRC6) ELISA Kit

CATALOG #
AAA137326

REACTIVITY
Human


PRODUCT SYNONYMS

Baculoviral IAP Repeat Containing Protein 6 (BIRC6); Human Baculoviral IAP Repeat Containing Protein 6 (BIRC6) ELISA Kit; BRUCE; Apollon; BIR repeat-containing ubiquitin-conjugating enzyme; Ubiquitin-conjugating BIR domain enzyme apollon; BIRC6 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Baculoviral IAP Repeat Containing Protein 6 (BIRC6). No significant cross-reactivity or interference between Baculoviral IAP Repeat Containing Protein 6 (BIRC6) and analogues was observed.
Sequence Length	4857
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.116ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Baculoviral IAP Repeat Containing Protein 6 (BIRC6) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Baculoviral IAP Repeat Containing Protein 6 (BIRC6) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BIRC6 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Baculoviral IAP Repeat Containing Protein 6 (BIRC6). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Baculoviral IAP Repeat Containing Protein 6 (BIRC6). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Baculoviral IAP Repeat Containing Protein 6 (BIRC6), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Baculoviral IAP Repeat Containing Protein 6 (BIRC6) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for BIRC6 elisa kit	Signal Transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	153792694
NCBI GeneID	57448
NCBI Accession #	NP_057336.3
NCBI GenBank Nucleotide #	NM_016252.3
UniProt Accession #	Q9NR09; Q9ULD1
Molecular Weight	530,255 Da
NCBI Official Full Name	baculoviral IAP repeat-containing protein 6
NCBI Official Synonym Full Names	baculoviral IAP repeat containing 6
NCBI Official Symbol	BIRC6
NCBI Official Synonym Symbols	BRUCE; APOLLON
NCBI Protein Information	baculoviral IAP repeat-containing protein 6
UniProt Protein Name	Baculoviral IAP repeat-containing protein 6
UniProt Gene Name	BIRC6
UniProt Synonym Gene Names	KIAA1289; BRUCE; APOLLON
UniProt Entry Name	BIRC6_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.