

Rat Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) ELISA Kit | ALPL elisa kit

Rat Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) ELISA Kit

CATALOG #
AAA137520

REACTIVITY
Rat


PRODUCT SYNONYMS

Alkaline Phosphatase; Liver/Bone/Kidney (ALPL); Rat Alkaline Phosphatase; Liver/Bone/Kidney (ALPL) ELISA Kit; BALP; BSAP; HOPS; AP-TNAP; TNSALP; Bone Alkaline Phosphatase; Tissue-Nonspecific Isozyme; ALPL elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	This assay has high sensitivity and excellent specificity for detection of Alkaline Phosphatase, Liver/Bone/Kidney (ALPL). No significant cross-reactivity or interference between Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) and analogues was observed.
Sequence Length	524
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	3.12-200ng/mL
Sensitivity	1.15ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ALPL elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Alkaline Phosphatase, Liver/Bone/Kidney (ALPL). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Alkaline Phosphatase, Liver/Bone/Kidney (ALPL). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Alkaline Phosphatase, Liver/Bone/Kidney (ALPL), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for ALPL elisa kit	Enzyme & Kinase; Metabolic pathway; Hepatology; Bone metabolism

NCBI AND UNIPROT INFORMATION

NCBI GI #	187954299
NCBI GeneID	249
Molecular Weight	51,045 Da
NCBI Official Full Name	Alkaline phosphatase, liver/bone/kidney
NCBI Official Synonym Full Names	alkaline phosphatase, liver/bone/kidney
NCBI Official Symbol	ALPL
NCBI Official Synonym Symbols	HOPS; TNAP; APTNAP; TNSALP; AP-TNAP
NCBI Protein Information	alkaline phosphatase, tissue-nonspecific isozyme
UniProt Protein Name	Alkaline phosphatase, tissue-nonspecific isozyme
UniProt Gene Name	ALPL
UniProt Synonym Gene Names	AP-TNAP; TNSALP
UniProt Entry Name	PPBT_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.