

Human Lysosomal Associated Membrane Protein 2 (LAMP2) ELISA Kit | LAMP2 elisa kit

Human Lysosomal Associated Membrane Protein 2 (LAMP2) ELISA Kit

CATALOG #
AAA138023

REACTIVITY
Human


PRODUCT SYNONYMS

Lysosomal Associated Membrane Protein 2 (LAMP2); Human Lysosomal Associated Membrane Protein 2 (LAMP2) ELISA Kit; CD107-B; CD107b; LAMP-2C; LAMPB; LGP96; Lysosome-associated membrane glycoprotein 2; CD107 antigen-like family member B; LAMP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Lysosomal Associated Membrane Protein 2 (LAMP2). No significant cross-reactivity or interference between Lysosomal Associated Membrane Protein 2 (LAMP2) and analogues was observed.
Sequence Length	365
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma and other biological fluids
Detection Range	62.5-4,000pg/mL
Sensitivity	24.9pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Lysosomal Associated Membrane Protein 2 (LAMP2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Lysosomal Associated Membrane Protein 2 (LAMP2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LAMP2 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Lysosomal Associated Membrane Protein 2 (LAMP2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Lysosomal Associated Membrane Protein 2 (LAMP2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Lysosomal Associated Membrane Protein 2 (LAMP2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Lysosomal Associated Membrane Protein 2 (LAMP2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for LAMP2 elisa kit	CD & Adhesion molecule; Infection immunity

NCBI AND UNIPROT INFORMATION

NCBI GI #	355698967
Molecular Weight	40,247 Da
NCBI Official Full Name	lysosomal-associated membrane protein 2, partial
UniProt Protein Name	Lysosomal-associated membrane protein 2
UniProt Entry Name	G9K7U8_MUSPF

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/lysosomal-associated-membrane-protein-2-lamp2/138023>