

Mouse Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit | TICAM1 elisa kit

Mouse Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit

CATALOG #
AAA138336

REACTIVITY
Mouse


PRODUCT SYNONYMS

Toll Like Receptor Adaptor Molecule 1 (TICAM1); Mouse Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit; TRIF; PRVTIRB; Toll-interleukin-1 receptor domain-containing adapter protein inducing interferon beta; Proline-rich; vinculin and TIR domain-containing protein B; TICAM1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Toll Like Receptor Adaptor Molecule 1 (TICAM1). No significant cross-reactivity or interference between Toll Like Receptor Adaptor Molecule 1 (TICAM1) and analogues was observed.
Sequence Length	712
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	62.5-4000pg/mL
Sensitivity	21.7pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Toll Like Receptor Adaptor Molecule 1 (TICAM1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Toll Like Receptor Adaptor Molecule 1 (TICAM1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TICAM1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Toll Like Receptor Adaptor Molecule 1 (TICAM1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Toll Like Receptor Adaptor Molecule 1 (TICAM1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Toll Like Receptor Adaptor Molecule 1 (TICAM1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Toll Like Receptor Adaptor Molecule 1 (TICAM1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for TICAM1 elisa kit	Immune molecule

NCBI AND UNIPROT INFORMATION

NCBI GI #	223460864
NCBI GeneID	148022
UniProt Accession #	Q8IUC6; O75532; Q86XP8; Q96GA0; B3Y691
Molecular Weight	76,422 Da
NCBI Official Full Name	Toll-like receptor adaptor molecule 1
NCBI Official Synonym Full Names	toll-like receptor adaptor molecule 1
NCBI Official Symbol	TICAM1
NCBI Official Synonym Symbols	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1
NCBI Protein Information	TIR domain-containing adapter molecule 1
UniProt Protein Name	TIR domain-containing adapter molecule 1
UniProt Gene Name	TICAM1
UniProt Synonym Gene Names	PRVTIRB; TRIF; TICAM-1; MyD88-3; TIR domain-containing adapter protein inducing IFN-beta
UniProt Entry Name	TCAM1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.