

Human Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) ELISA Kit | CC2D1A elisa kit

Human Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) ELISA Kit

CATALOG #
AAA138658

REACTIVITY
Human


PRODUCT SYNONYMS

Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A); Human Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) ELISA Kit; FREUD-1; MRT3; AKI1; Mental Retardation; Nonsyndromic; Autosomal Recessive 3; Five prime repressor element under dual repression-binding 1; Akt kinase-interacting protein 1; CC2D1A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A). No significant cross-reactivity or interference between Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) and analogues was observed.
Sequence Length	950
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.128ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CC2D1A elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for CC2D1A elisa kit	Signal Transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	96974985
NCBI GeneID	54862
NCBI Accession #	NP_060191.3
NCBI GenBank Nucleotide #	NM_017721.4
UniProt Accession #	Q6P1N0; Q7Z435; Q86XV0; Q8NF89; Q9H603; Q9NXI1
Molecular Weight	103,934 Da
NCBI Official Full Name	coiled-coil and C2 domain-containing protein 1A
NCBI Official Synonym Full Names	coiled-coil and C2 domain containing 1A
NCBI Official Symbol	CC2D1A
NCBI Official Synonym Symbols	MRT3; FREUD-1; Freud-1/Aki1
NCBI Protein Information	coiled-coil and C2 domain-containing protein 1A
UniProt Protein Name	Coiled-coil and C2 domain-containing protein 1A
UniProt Gene Name	CC2D1A
UniProt Synonym Gene Names	AKI1; FRE under dual repression-binding protein 1; Freud-1
UniProt Entry Name	C2D1A_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.