

Human Alkaline Phosphatase, Intestinal (ALPI) ELISA Kit | ALPI elisa kit

Human Alkaline Phosphatase, Intestinal (ALPI) ELISA Kit

CATALOG #
AAA139469

REACTIVITY
Human


PRODUCT SYNONYMS

Alkaline Phosphatase; Intestinal (ALPI); Human Alkaline Phosphatase; Intestinal (ALPI) ELISA Kit; CIP; IAP; Calf Intestinal Alkaline Phosphatase; ALPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Alkaline Phosphatase, Intestinal (ALPI). No significant cross-reactivity or interference between Alkaline Phosphatase, Intestinal (ALPI) and analogues was observed.
Sequence Length	528
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	1.56-100ng/mL
Sensitivity	0.63ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Alkaline Phosphatase, Intestinal (ALPI) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Alkaline Phosphatase, Intestinal (ALPI) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ALPI elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Alkaline Phosphatase, Intestinal (ALPI). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Alkaline Phosphatase, Intestinal (ALPI). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Alkaline Phosphatase, Intestinal (ALPI), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Alkaline Phosphatase, Intestinal (ALPI) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for ALPI elisa kit	Enzyme & Kinase; Metabolic pathway; Gastroenterology

NCBI AND UNIPROT INFORMATION

NCBI GI #	124376142
NCBI GeneID	248
Molecular Weight	56,812 Da
NCBI Official Full Name	Alkaline phosphatase, intestinal
NCBI Official Synonym Full Names	alkaline phosphatase, intestinal
NCBI Official Symbol	ALPI
NCBI Official Synonym Symbols	IAP
NCBI Protein Information	intestinal-type alkaline phosphatase
UniProt Protein Name	Intestinal-type alkaline phosphatase
UniProt Gene Name	ALPI
UniProt Synonym Gene Names	IAP; Intestinal alkaline phosphatase
UniProt Entry Name	PPBI_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.