

Human Prosaposin (PSAP) ELISA Kit | PSAP elisa kit

Human Prosaposin (PSAP) ELISA Kit

CATALOG #
AAA139591

REACTIVITY
Human


PRODUCT SYNONYMS

Prosaposin (PSAP); Human Prosaposin (PSAP) ELISA Kit; SAP1; GLBA; CSAct; Dispersin; Pro-Saposin; Proactivator polypeptide; Glucosylceramidase activator; Sulfatide/GM1 activator; Cerebroside sulfate activator; Sphingolipid Activator 1; PSAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Prosaposin (PSAP). No significant cross-reactivity or interference between Prosaposin (PSAP) and analogues was observed.
Sequence Length	527
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cerebrospinal fluid, Seminal plasma, Cell culture supernates and other biological fluids
Detection Range	0.78-50ng/mL
Sensitivity	0.33ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Prosaposin (PSAP) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Prosaposin (PSAP) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSAP elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Prosaposin (PSAP). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Prosaposin (PSAP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Prosaposin (PSAP), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Prosaposin (PSAP) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for PSAP elisa kit	Metabolic pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	11139234
NCBI GeneID	5660
UniProt Accession #	P07602; P07292; P15793; P78538; P78541; P78546; P78547; P78558; Q53Y86; Q6IBQ6; Q92739; Q92740
Molecular Weight	58,484 Da
NCBI Official Full Name	prosaposin, partial
NCBI Official Synonym Full Names	prosaposin
NCBI Official Symbol	PSAP
NCBI Official Synonym Symbols	GLBA; SAP1
NCBI Protein Information	prosaposin
UniProt Protein Name	Prosaposin
UniProt Gene Name	PSAP
UniProt Synonym Gene Names	GLBA; SAP1; CSAct; SAP-1; SAP-2
UniProt Entry Name	SAP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/prosaposin-psap/139591>