

Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit | BPI elisa kit

Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit

CATALOG #
AAA14009

REACTIVITY
Human


PRODUCT SYNONYMS

Bactericidal/Permeability Increasing Protein (BPI); Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit; CAP 57; BPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of BPI. No significant cross-reactivity or interference between BPI and analogues was observed.
Sequence Length	487
Assay Type	Quantitative Sandwich
Samples	Human serum, plasma, tissue homogenates or other biological fluids
Detection Range	12.5-800pg/mL
Sensitivity	< 5.4pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level BPI were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level BPI were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BPI elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to BPI. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to BPI. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain BPI, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of BPI in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of BPI in human serum, plasma, tissue homogenates or other biological fluids.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	157276599
NCBI GeneID	671
NCBI Accession #	NP_001716.2
NCBI GenBank Nucleotide #	NM_001725.2
UniProt Accession #	P17213; Q1ZZU8; Q5JRW0; Q8IW58; Q9BYZ9; Q9H1L2; Q9H1M8; Q9H203; Q9UCT4; Q9UD65; B2RCY2
Molecular Weight	53,900 Da
NCBI Official Full Name	bactericidal permeability-increasing protein
NCBI Official Synonym Full Names	bactericidal/permeability-increasing protein
NCBI Official Symbol	BPI
NCBI Official Synonym Symbols	rBPI; BPIFD1
NCBI Protein Information	bactericidal permeability-increasing protein
UniProt Protein Name	Bactericidal permeability-increasing protein
UniProt Gene Name	BPI
UniProt Synonym Gene Names	BPI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/bactericidal-permeability-increasing-protein-bpi/14009>