

## Rat Afamin (AFM) ELISA Kit | AFM elisa kit

Rat Afamin (AFM) ELISA Kit

CATALOG #  
AAA140110

REACTIVITY  
Rat


### PRODUCT SYNONYMS

Afamin (AFM); Rat Afamin (AFM) ELISA Kit; ALB2; ALBA; ALF; A-ALB; Alpha-Albumin; AFM elisa kit

### PRODUCT SPECIFICATIONS

|                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                                        | Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Specificity                                       | This assay has high sensitivity and excellent specificity for detection of Afamin (AFM). No significant cross-reactivity or interference between Afamin (AFM) and analogues was observed.                                                                                                                                                                                                                                                                                                  |
| Sequence Length                                   | 599                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Assay Type                                        | Competitive Inhibition                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Samples                                           | Serum, Plasma and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Detection Range                                   | 2.47-200ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Sensitivity                                       | 0.93ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Application                                       | Enzyme-linked immunosorbent assay for Antigen Detection.                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Intra-assay Precision (Precision within an assay) | 3 samples with low, middle and high level Afamin (AFM) were tested 20 times on one plate, respectively.                                                                                                                                                                                                                                                                                                                                                                                    |
| Inter-assay Precision (Precision between assays)  | 3 samples with low, middle and high level Afamin (AFM) were tested on 3 different plates, 8 replicates in each plate.                                                                                                                                                                                                                                                                                                                                                                      |
| CV(%) =                                           | SD/meanX100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Intra-Assay                                       | CV<10%                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Inter-Assay                                       | CV<12%                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Preparation and Storage                           | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for AFM elisa kit | The microtiter plate provided in this kit has been pre-coated with an antigen. Standards or samples are then added to the appropriate microtiter plate wells with a Horseradish Peroxidase (HRP)-conjugated secondary antibody. After TMB substrate solution is added, those wells that contain Afamin (AFM) will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Afamin (AFM) in the samples is then determined by comparing the O.D. of the samples to the standard curve. |
| Product Categories/Family for AFM elisa kit   | <a href="#">Tumor immunity</a> ; <a href="#">Infection immunity</a> ; <a href="#">Hepatology</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

NCBI AND UNIPROT INFORMATION

|                                  |                                |
|----------------------------------|--------------------------------|
| NCBI GI #                        | 4501987                        |
| NCBI GeneID                      | 173                            |
| NCBI Accession #                 | NP_001124.1                    |
| NCBI GenBank Nucleotide #        | NM_001133.2                    |
| UniProt Accession #              | P43652; Q32MR3; Q4W5C5; A8K3E1 |
| Molecular Weight                 | 69,069 Da                      |
| NCBI Official Full Name          | afamin                         |
| NCBI Official Synonym Full Names | afamin                         |
| NCBI Official Symbol             | AFM                            |
| NCBI Official Synonym Symbols    | ALF; ALB2; ALBA                |
| NCBI Protein Information         | afamin                         |
| UniProt Protein Name             | Afamin                         |
| UniProt Gene Name                | <b>AFM</b>                     |
| UniProt Synonym Gene Names       | ALB2; ALBA; Alpha-Alb          |
| UniProt Entry Name               | AFAM_HUMAN                     |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/afamin-afm/140110>