

Human Interferon Inducible Protein 35 (IFI35) ELISA Kit | IFI35 elisa kit

Human Interferon Inducible Protein 35 (IFI35) ELISA Kit

CATALOG #
AAA140691

REACTIVITY
Human


PRODUCT SYNONYMS

Interferon Inducible Protein 35 (IFI35); Human Interferon Inducible Protein 35 (IFI35) ELISA Kit; IFP35; IFI35 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Interferon Inducible Protein 35 (IFI35). No significant cross-reactivity or interference between Interferon Inducible Protein 35 (IFI35) and analogues was observed.
Sequence Length	288
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	78-5,000pg/mL
Sensitivity	33pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Interferon Inducible Protein 35 (IFI35) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Interferon Inducible Protein 35 (IFI35) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IFI35 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Interferon Inducible Protein 35 (IFI35). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Interferon Inducible Protein 35 (IFI35). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Interferon Inducible Protein 35 (IFI35), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Interferon Inducible Protein 35 (IFI35) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for IFI35 elisa kit	Cytokine

NCBI AND UNIPROT INFORMATION

NCBI GI #	306922358
NCBI GeneID	3430
NCBI Accession #	NP_005524.2
NCBI GenBank Nucleotide #	NM_005533.4
Molecular Weight	31,777 Da
NCBI Official Full Name	interferon-induced 35 kDa protein
NCBI Official Synonym Full Names	interferon-induced protein 35
NCBI Official Symbol	IFI35
NCBI Official Synonym Symbols	IFP35
NCBI Protein Information	interferon-induced 35 kDa protein
UniProt Protein Name	Interferon-induced 35 kDa protein
UniProt Gene Name	IFI35
UniProt Synonym Gene Names	IFP35; IFP 35; Ifi-35
UniProt Entry Name	IN35_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/interferon-inducible-protein-35-ifi35/140691>