

Gamma-Aminobutyric Acid (gABA) DIY ELISA Kit | gABA diy elisa kit

General Gamma-Aminobutyric Acid (gABA) ELISA Kit DIY Materials

CATALOG # AAA147665	REACTIVITY Pan-species (General)
-------------------------------	--


PRODUCT SYNONYMS

Gamma-Aminobutyric Acid (gABA); General Gamma-Aminobutyric Acid (gABA) ELISA Kit DIY Materials; 4-Aminobutyric Acid; gABA diy elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Pan-species (General)
Specificity	This assay has high sensitivity and excellent specificity for detection of E2.
Sequence Length	462
Assay Type	Competitive
Detection Range	12.35-1,000pg/mL.
Preparation and Storage	Antibodies, Standard and Streptavidin-HRP should be stored at -20 degree C. TMB should be stored at 4 degree C. 96-well Plate could be stored at room temperature. The contents are valid for twelve months. They are stable for one month after opening when stored at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for gABA diy elisa kit	Intended Uses: This 'Do it Yourself (DIY)' assay kit contains materials for developing your own immunoassays. This may include the development of competitive inhibition ELISA to measure E2 in serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids in vitro. This kit contains sufficient materials for preparation of at least ten 96-well plates. Sample Preparation, Assay Protocol, and Detection Range for Reference are provided as suggestions only. Researchers should optimize the use of kit materials and protocols for their own model system and determine if the kit is suitable for their samples and immunoassay.
Product Categories/Family for gABA diy elisa kit	Signal transduction ; Cardiovascular biology ; Neuroscience ; Bone metabolism

NCBI AND UNIPROT INFORMATION

NCBI GI #	60654715
NCBI Official Full Name	gamma-aminobutyric acid

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/general-diy-elisa-kits/gamma-aminobutyric-acid-gaba/147665>