

25-Hydroxyvitamin D3 (HVD3) DIY ELISA Kit | HVD3 diy elisa kit

General 25-Hydroxyvitamin D3 (HVD3) ELISA Kit DIY Materials

CATALOG #
AAA147668

REACTIVITY
Pan-species (General)


PRODUCT SYNONYMS

25-Hydroxyvitamin D3 (HVD3); General 25-Hydroxyvitamin D3 (HVD3) ELISA Kit DIY Materials; 25-HVD3; 25HVD3; Calcifediol; Calcidiol; H-VD3; 25-Hydroxycholecalciferol; 25-Hydroxy Vitamin D3; 25(OH)D; HVD3 diy elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Pan-species (General)
Specificity	This assay has high sensitivity and excellent specificity for detection of HVD3.
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Competitive
Detection Range	246.9-20,000pg/mL
Preparation and Storage	Antibodies, Standard and Streptavidin-HRP should be stored at -20 degree C. TMB should be stored at 4 degree C. 96-well Plate could be stored at room temperature. The contents are valid for twelve months. They are stable for one month after opening when stored at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HVD3 diy elisa kit

Intended Uses: This 'Do it Yourself (DIY)' assay kit contains materials for developing your own immunoassays. This may include the development of competitive inhibition ELISA to measure HVD3 in serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids in vitro. This kit contains sufficient materials for preparation of at least ten 96-well plates. Sample Preparation, Assay Protocol, and Detection Range for Reference are provided as suggestions only. Researchers should optimize the use of kit materials and protocols for their own model system and determine if the kit is suitable for their samples and immunoassay.

Product Categories/Family for HVD3 diy elisa kit

[Reproductive science](#); [Genetic science](#); [Nutrition metabolism](#); [Bone metabolism](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.