

Rat Tight Junction Associated Protein 1, Peripheral (TJAP1) ELISA Kit | TJAP1 elisa kit

Rat Tight Junction Associated Protein 1, Peripheral (TJAP1) ELISA Kit

CATALOG #
AAA148653

REACTIVITY
Rat


PRODUCT SYNONYMS

Tight Junction Associated Protein 1; Peripheral (TJAP1); Rat Tight Junction Associated Protein 1; Peripheral (TJAP1) ELISA Kit; PILT; TJP4; Tight Junction Protein 4; Protein incorporated later into tight junctions; TJAP1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	This assay has high sensitivity and excellent specificity for detection of Tight Junction Associated Protein 1, Peripheral (TJAP1). No significant cross-reactivity or interference between Tight Junction Associated Protein 1, Peripheral (TJAP1) and analogues was observed.
Sequence Length	547
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	78-5,000pg/mL
Sensitivity	34pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Tight Junction Associated Protein 1, Peripheral (TJAP1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Tight Junction Associated Protein 1, Peripheral (TJAP1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TJAP1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Tight Junction Associated Protein 1, Peripheral (TJAP1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Tight Junction Associated Protein 1, Peripheral (TJAP1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Tight Junction Associated Protein 1, Peripheral (TJAP1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Tight Junction Associated Protein 1, Peripheral (TJAP1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for TJAP1 elisa kit	Signal Transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	119624601
NCBI GeneID	93643
Molecular Weight	53,174 Da
NCBI Official Full Name	tight junction associated protein 1 (peripheral), isoform CRA_a
NCBI Official Synonym Full Names	tight junction associated protein 1
NCBI Official Symbol	TJAP1
NCBI Official Synonym Symbols	PILT; TJP4
NCBI Protein Information	tight junction-associated protein 1
UniProt Protein Name	Tight junction-associated protein 1
UniProt Gene Name	TJAP1
UniProt Synonym Gene Names	PILT; TJP4

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/tight-junction-associated-protein-1-peripheral-tjap1/148653>