

Human Phosphoglycerate Mutase 5 (PGAM5) ELISA Kit | PGAM5 elisa kit

Human Phosphoglycerate Mutase 5 (PGAM5) ELISA Kit

CATALOG #
AAA148699

REACTIVITY
Human


PRODUCT SYNONYMS

Phosphoglycerate Mutase 5 (PGAM5); Human Phosphoglycerate Mutase 5 (PGAM5) ELISA Kit; BXLbv68; Bcl-XL-binding protein v68; Serine/threonine-protein phosphatase PGAM5; mitochondrial; PGAM5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Phosphoglycerate Mutase 5 (PGAM5). No significant cross-reactivity or interference between Phosphoglycerate Mutase 5 (PGAM5) and analogues was observed.
Sequence Length	288
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue Homogenates And Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	0.057ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Phosphoglycerate Mutase 5 (PGAM5) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Phosphoglycerate Mutase 5 (PGAM5) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGAM5 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Phosphoglycerate Mutase 5 (PGAM5). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Phosphoglycerate Mutase 5 (PGAM5). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Phosphoglycerate Mutase 5 (PGAM5), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Phosphoglycerate Mutase 5 (PGAM5) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for PGAM5 elisa kit	Enzyme & Kinase

NCBI AND UNIPROT INFORMATION

NCBI GI #	665389201
NCBI GeneID	31143
NCBI Accession #	NP_001284795.1
NCBI GenBank Nucleotide #	NM_001297866.1
Molecular Weight	33,166 Da
NCBI Official Full Name	phosphoglycerate mutase 5, isoform B
NCBI Official Synonym Full Names	Phosphoglycerate mutase 5
NCBI Official Symbol	Pgam5
NCBI Official Synonym Symbols	CG14816; DmelCG14816; dPGAM5; DPGAM5; EG:63B12.4; PGAM5
NCBI Protein Information	CG14816 gene product from transcript CG14816-RB
UniProt Protein Name	Serine/threonine-protein phosphatase Pgam5, mitochondrial
UniProt Gene Name	Pgam5
UniProt Synonym Gene Names	DPGAM5

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/phosphoglycerate-mutase-5-pgam5/148699>