

Human Elastase 2, Neutrophil (ELA2) Antibody Pair Kit | ELA2 antibody pair kit

Human Elastase 2, Neutrophil (ELA2) Antibody Pair Kit (with Standard)

CATALOG # AAA148729	REACTIVITY Human	APPLICATIONS ELISA, ELISA
-------------------------------	----------------------------	-------------------------------------


PRODUCT SYNONYMS

Elastase 2; Neutrophil (ELA2); Human Elastase 2; Neutrophil (ELA2) Antibody Pair Kit (with Standard); ELANE; HLE; HNE; NE; PMN-E; Neutrophil Elastase; Medullasin; Polymorphonuclear Leukocyte Elastase; Bone Marrow Serine Protease; Human Leukocyte Elastase; PMN Elastase; ELA2 antibody pair kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, Plasma, Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.38-24ng/mL
Applicable Applications for ELA2 antibody pair kit	Luminex, ELISPOT, CLIA (Chemiluminescence Immunoassay), ELISA
Standard	Recombinant ELA2 (Lyophilized, 1 vial-10ug)
Capture Antibody	Rabbit polyclonal antibody (Liquid, 1 vial-400μg/800μL)
Detection Antibody	Biotinylated Rabbit polyclonal antibody (Liquid, 1 vial-100μg/500μL)
Product Use	Raw materials can be used for preparation of 5x96T or 10x96T ELISA Kits.
Antibody Support Pack	Antibody Support Pack (Buffer Pack) is sold separately as

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Product Categories/Family for ELA2 antibody pair kit	Enzyme & Kinase; Infection immunity
--	-------------------------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	46361516
NCBI GeneID	1991
Molecular Weight	28,518 Da
NCBI Official Full Name	elastase 2, neutrophil

NCBI Official Synonym Full Names	elastase, neutrophil expressed
NCBI Official Symbol	ELANE
NCBI Official Synonym Symbols	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E
NCBI Protein Information	neutrophil elastase
UniProt Protein Name	Neutrophil elastase
UniProt Gene Name	ELANE
UniProt Synonym Gene Names	ELA2; HLE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-antibody-pair-kits/elastase-2-neutrophil-ela2/148729>