

Bovine Reactive Oxygen Species Modulator 1 (ROMO1) ELISA Kit | ROMO1 elisa kit

Bovine Reactive Oxygen Species Modulator 1 (ROMO1) ELISA Kit

CATALOG #
AAA149209

REACTIVITY
Bovine


PRODUCT SYNONYMS

Reactive Oxygen Species Modulator 1 (ROMO1); Bovine Reactive Oxygen Species Modulator 1 (ROMO1) ELISA Kit; C20orf52; MTGMP; Glyrichin; Mitochondrial Targeting GXXXG Protein; Epididymis tissue protein Li 175; Protein MGR2 homolog; ROMO1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Specificity	This assay has high sensitivity and excellent specificity for detection of ROBO1. No significant cross-reactivity or interference between ROBO1 and analogues was observed.
Sequence Length	79
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates And Other Biological Fluids
Detection Range	0.312-20ng/mL
Sensitivity	0.119ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level ROBO1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level ROBO1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ROMO1 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of ROBO1 in human tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to ROBO1. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to ROBO1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain ROBO1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of ROBO1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for ROMO1 elisa kit	Signal Transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	18152785
NCBI GeneID	140823
NCBI Accession #	NP_542786.1
NCBI GenBank Nucleotide #	NM_080748.3
UniProt Accession #	P60602; Q3MHD5; Q5QP16; Q9CQ98; Q9H1N2; A7M872; E1P5R9; E9KL28
Molecular Weight	5,900 Da
NCBI Official Full Name	reactive oxygen species modulator 1
NCBI Official Synonym Full Names	reactive oxygen species modulator 1
NCBI Official Symbol	ROMO1
NCBI Official Synonym Symbols	MTGM; MTGMP; C20orf52; bA353C18.2
NCBI Protein Information	reactive oxygen species modulator 1
UniProt Protein Name	Reactive oxygen species modulator 1
UniProt Gene Name	ROMO1
UniProt Synonym Gene Names	C20orf52; ROS modulator 1; MTGM

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/reactive-oxygen-species-modulator-1-romo1/149209>