

Human Cartilage Intermediate Layer Protein 2 (CILP2) CLIA Kit | CILP2 clia kit

Human Cartilage Intermediate Layer Protein 2 (CILP2)

CATALOG #
AAA149988

REACTIVITY
Human


PRODUCT SYNONYMS

Cartilage Intermediate Layer Protein 2 (CILP2); Human Cartilage Intermediate Layer Protein 2 (CILP2); Cartilage Intermediate Layer Protein 2; CILP2 clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Cartilage Intermediate Layer Protein 2 (CILP2). No significant cross-reactivity or interference between Cartilage Intermediate Layer Protein 2 (CILP2) and analogues was observed.
Sequence Length	1156
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates and other biological fluids
Detection Range	27.4-20,000pg/mL
Sensitivity	10.1pg/mL
Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Cartilage Intermediate Layer Protein 2 (CILP2) were tested 20 times on one plate, respectively. Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Cartilage Intermediate Layer Protein 2 (CILP2) were tested on 3 different plates, 8 replicates in each plate. $CV(\%) = SD/mean \times 100$ Intra-Assay: $CV < 10\%$ Inter-Assay: $CV < 12\%$
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CILP2 clia kit	The microplate provided in this kit has been pre-coated with an antibody specific to Cartilage Intermediate Layer Protein 2 (CILP2). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Cartilage Intermediate Layer Protein 2 (CILP2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Cartilage Intermediate Layer Protein 2 (CILP2) level in the sample or standard.
Product Categories/Family for CILP2 clia kit	Metabolic Pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	145699139
NCBI Accession #	NP_694953.2
NCBI GenBank Nucleotide #	NP_694953.2
UniProt Accession #	Q8IUL8
NCBI Official Full Name	cartilage intermediate layer protein 2
UniProt Protein Name	Cartilage intermediate layer protein 2
UniProt Gene Name	CILP2
UniProt Synonym Gene Names	CILP-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/cartilage-intermediate-layer-protein-2-cilp2/149988>