

Human Cartilage Intermediate Layer Protein (CILP) ELISA Kit | CILP elisa kit

Human Cartilage Intermediate Layer Protein (CILP) ELISA Kit

CATALOG #
AAA150341

REACTIVITY
Human


PRODUCT SYNONYMS

Cartilage Intermediate Layer Protein (CILP); Human Cartilage Intermediate Layer Protein (CILP) ELISA Kit; CILP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of CILP. No significant cross-reactivity or interference between CILP and analogues was observed.
Samples	Human tissue homogenates, cell lysates and other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	< 0.054ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level CILP were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level CILP were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CILP elisa kit

Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of CILP in human tissue homogenates, cell lysates and other biological fluids.

Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to CILP. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to CILP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain CILP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of CILP in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cartilage-intermediate-layer-protein-cilp/150341>