

## Human Lysine Specific Demethylase 3A (KDM3A) ELISA Kit | KDM3A elisa kit

Human Lysine Specific Demethylase 3A (KDM3A) ELISA Kit

|                               |                            |                              |
|-------------------------------|----------------------------|------------------------------|
| CATALOG #<br><b>AAA150526</b> | REACTIVITY<br><b>Human</b> | APPLICATIONS<br><b>ELISA</b> |
|-------------------------------|----------------------------|------------------------------|


### PRODUCT SYNONYMS

Lysine Specific Demethylase 3A (KDM3A); Human Lysine Specific Demethylase 3A (KDM3A) ELISA Kit; JMJD1; JMJD1A; TSGA; JHMD2A; JmjC domain-containing histone demethylation protein 2A; Jumonji domain-containing protein 1A; KDM3A elisa kit

### PRODUCT SPECIFICATIONS

|                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                                  | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity                                 | This assay has high sensitivity and excellent specificity for detection of KDM3A. No significant cross-reactivity or interference between KDM3A and analogues was observed.                                                                                                                                                                                                                                                                                                                |
| Assay Type                                  | Double-antibody Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Samples                                     | Tissue Homogenates, Cell Lysates And Other Biological Fluids                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Detection Range                             | 0.156-10ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Sensitivity                                 | 0.057ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Applicable Applications for KDM3A elisa kit | ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Intra-assay Precision                       | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level KDM3A were tested 20 times on one plate, respectively. Intra-Assay: CV<10%                                                                                                                                                                                                                                                                                                                    |
| Inter-assay Precision                       | Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level KDM3A were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%                                                                                                                                                                                                                                                                                  |
| Preparation and Storage                     | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for KDM3A elisa kit | <p>Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of KDM3A in human tissue homogenates, cell lysates and other biological fluids.</p> <p>Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to KDM3A. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to KDM3A. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain KDM3A, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of KDM3A in the samples is then determined by comparing the O.D. of the samples to the standard curve.</p> |
| Product Categories/Family for KDM3A elisa kit   | Enzyme & Kinase                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/lysine-specific-demethylase-3a-kdm3a/150526>