

Human Nuclear matrix protein 22, NMP-22 ELISA Kit | NMP-22 elisa kit

Human Nuclear matrix protein 22, NMP-22 ELISA Kit

CATALOG #
AAA15130

REACTIVITY
Human


PRODUCT SYNONYMS

Nuclear matrix protein 22; NMP-22; Human Nuclear matrix protein 22; NMP-22 ELISA Kit; NMP-22 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human NMP-22. No significant cross-reactivity or interference between human NMP-22 and analogues was observed.
Assay Type	Competitive
Samples	Serum, urine
Detection Range	0.156 ng/ml-40 ng/ml
Sensitivity	< 0.156ng/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<12%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<15%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NMP-22 elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with an antigen. Standards or samples are added to the appropriate microtiter plate wells with antibody specific for NMP-22 and Horseradish Peroxidase (HRP) conjugated goat-anti-mouse antibody. The competitive inhibition reaction is launched between with pre-coated NMP-22 and NMP-22 in samples. A substrate solution is added to the wells and the color develops in opposite to the amount of NMP-22 in the sample. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/nuclear-matrix-protein-22-nmp-22/15130>