

Rat vasoactive intestinal peptide ELISA Kit | VIP elisa kit

Rat Vasoactive Intestinal Peptide, VIP ELISA Kit

CATALOG #
AAA15209

REACTIVITY
Rat


PRODUCT SYNONYMS

vasoactive intestinal peptide; Rat Vasoactive Intestinal Peptide; VIP ELISA Kit; MGC13587; PHM27; ; VIP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	This assay has high sensitivity and excellent specificity for detection of rat VIP. No significant cross-reactivity or interference between rat VIP and analogues was observed.
Sequence Length	72
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and tissue homogenates
Detection Range	1.56 pg/ml-100 pg/ml
Sensitivity	< 0.39 pg/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VIP elisa kit	Principle of the Assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for VIP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VIP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VIP is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VIP bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	418255
UniProt Accession #	P32649
Molecular Weight	8,198 Da
NCBI Official Full Name	VIP peptides
UniProt Protein Name	VIP peptides
UniProt Gene Name	VIP
UniProt Synonym Gene Names	VIP
UniProt Entry Name	VIP_RABIT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/vip-rat-elisa-kits/vasoactive-intestinal-peptide/15209>